


Product Data Sheet

Cas No.:	1246819-06-6	Cat. No:	PC70945
Product Name:	兰索拉唑杂质10对照品		
Product synonym:	1-甲基-12-硫氧代吡啶并[1,2,3,4]咪唑-[1,2-a]苯并咪唑-2(12H)-酮;2-氧代戊烷-1,5-二甲酸二乙酯;兰索拉唑硫代杂质;硫化钬;兰索拉唑杂质10		
Chemical name:	兰索拉唑杂质10对照品		
MF:	C14H9N3OS	FW:	267.30576
Purity:	≥95%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C1C=CN2C3=NC4=CC=CC=C4N3C(=S)C2=C1C		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

仅用于研发试验中鉴别、检查。1246819-06-6|兰索拉唑杂质10对照品，制备中，详情请联系Medlife客服进行咨询。1-methyl-12-thioxobenzo[4',5']imidazo[2',1':2,3]imidazo[1,5-a]pyridin-2(12H)-one可放大生产，批量保障。