


Product Data Sheet

Cas No.:	116368-90-2	Cat. No:	PC58158
Product Name:	Acetylexidonin		
Product synonym:	-		
Chemical name:	Acetylexidonin		
MF:	C26H34O9	FW:	490.5428
Purity:	98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O1C(C23C(C(=C([H])[H])C([H])(C([H])([H])C([H])(C2([H])C2(C1([H])[H])C([H])(C([H])([H])C([H])([H])C(C([H])([H])[H])(C([H])([H])[H])C2([H])C([H])([H])OC(C([H])([H])[H])=O)OC(C([H])([H])[H])=O)OC(C([H])([H])[H])=O)C3([H])[H])=O)=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Acetylexidonin 是一种具有抗炎和抗癌活性的二萜化合物。Acetylexidonin 抑制肿瘤细胞的 IC₅₀ 分别为 3.69 μM (NB4) 和 26.22 μM (SHSY5Y)。