


Product Data Sheet

Cas No.:	549-21-3	Cat. No:	PC58176
Product Name:	Oxyberberine		
Product synonym:	氧基白白林;氧化小檗碱;小檗浸碱;氧化小檗碱;8-氧化小檗碱		
Chemical name:	Oxyberberine		
MF:	C20H17NO5	FW:	351.3527
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O1C([H])([H])OC2=C1C([H])=C1C(=C2[H])C2=C([H])C3C([H])=C([H])C(=C(C=3N2C([H])([H])C1([H])([H])=O)OC([H])([H])([H])OC([H])([H])([H])[H])		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Oxyberberine (Oxyberberin) 是一种天然生物碱类物质。Oxyberberine (Oxyberberin) 是一种天然生物碱类物质。

生物活性	Oxyberberine (Oxyberberin) is a natural alkaloid isolated from many plants.
包装储存	4°C, Keep In Dark; *In solvent : -80°C, 6 months; -20°C, 1 month (Keep In Dark);
溶解度数据	体外研究: DMSO : 16.67 mg/mL(47.45 mM);ultrasonic and warming and heat to 60°C) 配制储存液