


Product Data Sheet

Cas No.:	4375-07-9	Cat. No:	PC57979
Product Name:	(-)-Epipodophyllotoxin		
Product synonym:	-		
Chemical name:	(-)-Epipodophyllotoxin		
MF:	C22H22O8	FW:	414.40528
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	COC1C(OC)=CC([C@@H]2C3=CC4OCOC=4C=C3[C@@H](O)[C@H]3COC(=O)[C@H]23)=CC=1OC		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

(-)-Epipodophyllotoxin 是一种从 *Podophyllum peltatum* 中分离的癌细胞抗增殖剂，在 HeLa 细胞和 MCF-7 细胞中的 GI50 值分别为 0.36 和 0.24 μM。(-)-Epipodophyllotoxin 可在体外抑制有丝分裂纺锤体组装。