


Product Data Sheet

Cas No.:	486-28-2	Cat. No:	PC57976
Product Name:	Fraxinol		
Product synonym:	白蜡树精;白蜡树精,fraxinol,植物提取物,标准品,对照品;白蜡树精对照品;白藜芦醇;皮酚;秦皮醇		
Chemical name:	Fraxinol		
MF:	C11H10O5	FW:	222.1941
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O1C(C([H])=C([H])C2=C(C(=C(C([H])=C12)OC([H])([H])([H])O[H])OC([H])([H])([H])=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Fraxinol 是从半边莲 *Lobelia chinensis* 中分离出来的。Fraxinol 是从半边莲 *Lobelia chinensis* 中分离出来的。

生物活性	Fraxinol is isolated from <i>Lobelia chinensis</i> .
包装储存	4°C, Sealed and stored away from moisture and light; *In solvent : -80°C, 6 months; -20°C, 1 month (Sealed and stored away from moisture and light);
溶解度数据	体外研究: DMSO : 25 mg/mL(112.52 mM);ultrasonic and warming and heat to 60°C) 配制储存液