


Product Data Sheet

Cas No.:	157478-01-8	Cat. No:	PC63226
Product Name:	Yubeinine		
Product synonym:	裕贝甲素		
Chemical name:	Yubeinine		
MF:	C27H43NO3	FW:	429.63500
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	[C@@]12([H])CC[C@H](C)CN1C[C@@]1([C@](CC[C@]3([H])[C@@]1([H])C[C@@]1([H])[C@@]3([H])CC(=O)[C@]3([H])[C@@]1(C[C@@H](O)C3C)([C@@]2(O)C)[H])[H]		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Yubeinine 是一种生物碱，具有气管松弛作用。Yubeinine 是一种生物碱，具有气管松弛作用。

生物活性	Yubeinine is an alkaloid with tracheal relaxant effects.
包装储存	4°C, Keep In Dark; *In solvent : -80°C, 6 months; -20°C, 1 month (Keep In Dark);
溶解度数据	体外研究: DMSO : 100 mg/mL(232.75 mM;Need ultrasonic) 配制储存液