


Product Data Sheet

Cas No.:	1075281-70-7	Cat. No:	PL01730
Product Name:	Teslexivir hydrochloride		
Product synonym:	-		
Chemical name:	Teslexivir hydrochloride		
MF:	C35H37BRCLN3O4	FW:	679.042987585068
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	N(C1C=CC(CC2C=CC=CC=2)=CC=1)(C1C(=CC=CC=1N1CCCCC1)C)NC(=O)CC1C(=CC(C(=O)O)=C(OC)C=1)Br.Cl		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Teslexivir (BTA074) hydrochloride 是一种有效的抗病毒剂。Teslexivir hydrochloride 是两种关键病毒蛋白 (E1 和 E2) 相互作用的强效选择性抑制剂, 这两种蛋白的相互作用是人乳头瘤病毒 (HPV) 6 和 11 的 DNA复制和病毒产生的重要步骤。Teslexivir hydrochloride 可用于尖锐湿疣的研究。

生物活性	Teslexivir (BTA074) hydrochloride is a potent antiviral agent. Teslexivir hydrochloride is a potent and selective inhibitor of the interaction between two essential viral proteins, E1 and E2, an association that is a necessary step in the DNA replication and thus viral production for Human Papilloma Virus (HPV) 6 and 11. Teslexivir hydrochloride can be used for condyloma research.
包装储存	4°C, sealed storage, away from moistur In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)
溶解度数据	In Vitro: DMSO : 250 mg/mL (368.17 mM; Need ultrasonic)配制储备液