


Product Data Sheet

Cas No.:	2561476-24-0	Cat. No:	PL03738
Product Name:	EMAC10101d		
Product synonym:	-		
Chemical name:	EMAC10101d		
MF:	C17H15CL2N3O2S2	FW:	428.355899095535
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	ClC1C=C(C=CC=1C1=CS/C(=N\C2C=CC(=CC=2)S(N)(=O)=O)/N1CC)Cl		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

EMAC10101d 是有效的、选择性的 hCA II 抑制剂，其 K_i 值为 8.1 nM。

生物活性	EMAC10101d is a potent and selective toward hCA II inhibitor, with a K _i of 8.1 nM.
体外研究(In Vitro)	EMAC10101d exhibits K _i values of 9627.4 nM, 8.1 nM, 224.6 nM and 154.9 nM for hCA I , hCA II, hCA IX and hCA XII, respectively. has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 62.5 mg/mL (145.91 mM; Need ultrasonic)配制储备液