


Product Data Sheet

Cas No.:	1187568-16-6	Cat. No:	PL03995
Product Name:	DK3		
Product synonym:	7-methoxy-3-methyl-2-phenylquinazolin-4-one		
Chemical name:	DK3		
MF:	C16H14N2O2	FW:	266.294563770294
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	N1C2=C(C=CC(OC)=C2)C(=O)N(C)C=1C1=CC=CC=C1		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

DK3 是一种有效选择性的雌激素相关受体 α (ERRα) 激动剂。

生物活性	DK3 is a potent and selective estrogen-related receptor alpha (ERRα) agonist.
包装储存	Powder -20°C 3 years; In solvent -80°C 6 months
溶解度数据	In Vitro: DMSO : 12.5 mg/mL (46.94 mM; ultrasonic and warming and heat to 60°C)配制储备液