


Product Data Sheet

Cas No.:	1253584-63-2	Cat. No:	PL04259
Product Name:	Hsp90-IN-17 hydrochloride		
Product synonym:	-		
Chemical name:	Hsp90-IN-17 hydrochloride		
MF:	C21H21CLN4O7	FW:	476.867044210434
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O(C1C=CC(N(=O)=O)=CC=1)C1C=C(O)C=C(O)C=1C1ON=C(C(=O)NC2CCNCC2)C=1.Cl		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Hsp90-IN-17 (Example 5) hydrochloride 是一种 HSP90 抑制剂，可用于增殖性疾病，例如癌症和神经退行性疾病的研究。

生物活性	Hsp90-IN-17 (Example 5) hydrochloride is an HSP90 inhibitor that can be used in the study of proliferative diseases, such as cancer and neurodegenerative diseases.
包装储存	4°C, sealed storage, away from moisture and light In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)
溶解度数据	In Vitro: DMSO : 100 mg/mL (209.70 mM; Need ultrasonic)配制储备液