


Product Data Sheet

Cas No.:	2471-70-7	Cat. No:	PL05570
Product Name:	6-Methoxy-2-naphthoic acid		
Product synonym:	6-甲氧基-2-萘甲酸;6-甲氧萘-2羧酸;2-羧基-6-甲氧基萘 标准品;6-Methoxy-2-naphthoic Acid 6-甲氧基-2-萘甲酸;6-甲氧基萘-2-羧酸;萘普生EP杂质;萘普生相关化合物A USP标准品;萘普生杂质A;萘普生EP杂质(2-羧基-6-甲氧基萘)		
Chemical name:	6-Methoxy-2-naphthoic acid		
MF:	C12H10O3	FW:	202.2060
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C([H])([H])[H]C1C([H])=C([H])C2C([H])=C(C(=O)O[H])C([H])=C([H])C=2C=1[H]		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

6-Methoxy-2-naphthoic acid 是一个 NMDA 受体调节剂来自专利 WO 2012019106 A2。

生物活性	6-Methoxy-2-naphthoic acid is an NMDA receptor modulator extracted from patent WO 2012019106 A2.
IC ₅₀ & Target[1][2]	NMDA receptor
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 60 mg/mL (296.72 mM; Need ultrasonic and warming)配制储备液