

Product Data Sheet

Cas No.:		Cat. No:	PL05876
Product Name:	LCMV gp33-41 TFA		
Product synonym:	-		
Chemical name:	LCMV gp33-41 TFA		
MF:	-	FW:	-
Purity:	≥99%	Batch No.:	-
Storage:			
λmax:	-	Formulation:	-
Solubility :	-		
SMILES :	-		
InChI Code:	-		
InChI Key:	-		
WARNING This product is not for human or veterinary use.			

Product Description

LCMV gp33-41 (TFA), 羧基延伸的 11-aa-长肽, 受 MHC I 类 H-2Db 分子限制的淋巴细胞脉络丛脑膜炎病毒序列, 可呈递给细胞毒性 T 淋巴细胞。

生物活性	LCMV gp33-41 (TFA), the carboxyl-extended 11-aa-long peptide, is an lymphocytic choriomeningitis virus sequence restricted by MHC class I H-2Db molecules and presented to cytotoxic T lymphocytes.
包装储存	Sealed storage, away from moisture Powder -80°C 2 years; -20°C 1 year
溶解度数据	In Vitro: DMSO : 100 mg/mL (86.34 mM; Need ultrasonic)配制储备液