


Product Data Sheet

Cas No.:	667899-73-2	Cat. No:	PL05924
Product Name:	Histone H4 (2-21)		
Product synonym:	-		
Chemical name:	Histone H4 (2-21)		
MF:	C82H150N36O22	FW:	1992.29641389847
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	C(O)(=O)[C@H](CCCCN)NC(=O)[C@H](CCCNC(N)=N)NC(=O)[C@H](CC1N=CNC=1)NC(=O)[C@H](CCCNC(N)=N)NC(=O)[C@H](CCCCN)NC(=O)[C@H](C)NC(=O)CNC(=O)CNC(=O)[C@H](CCCCN)NC(=O)CNC(=O)[C@H](CC(C)C)NC(=O)CNC(=O)[C@H](CCCCN)NC(=O)CNC(=O)CNC(=O)[C@H](CCCCN)NC(=O)CNC(=O)[C@H](CCCNC(N)=N)NC(=O)CNC(=O)[C@H](CO)N		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Histone H4 (2-21) 是一组与单纯疱疹病毒1 (HSV-1) 基因组染色化相关的核心组蛋白。

生物活性	Histone H4 (2-21) is the core histones associated with chromatinization of herpes simplex virus 1 (HSV-1) genomes.
IC50 & Target[1][2]	HSV-1
包装储存	Sealed storage, away from moisturePowder -80°C 2 years; -20°C 1 year
溶解度数据	In Vitro: H ₂ O : 100 mg/mL (50.19 mM; Need ultrasonic)配制储备液