


Product Data Sheet

Cas No.:	142356-33-0	Cat. No:	PL06336
Product Name:	Boc-NH-C6-Br		
Product synonym:	N-BOC-6-溴代己胺;6-(叔丁氧羰基氨基)溴己烷;N-BOC-6-溴-己胺;(6-溴己基)氨基甲酸叔丁酯;6-(Boc-氨基)溴己烷;N-(6-溴己基)氨基甲酸叔丁酯		
Chemical name:	Boc-NH-C6-Br		
MF:	C11H22BRNO2	FW:	280.2019
Purity:	≥95%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	BrC([H])([H])C([H])([H])C([H])([H])C([H])([H])C([H])([H])C([H])([H])N([H])C(=O)OC(C([H])([H])[H])(C([H])([H])[H])C([H])([H])[H]		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Boc-NH-C6-Br 是一种不可降解 (non-cleavable) 的 linker，用于抗体偶联药物 (ADC)。

生物活性	Boc-NH-C6-Br is a non-cleavable linker used for antibody-drug conjugates (ADC).
IC50 & Target[1][2]	Non-cleavable
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 40 mg/mL (142.76 mM; Need ultrasonic)H ₂ O : 1 mg/mL (3.57 mM; ultrasonic and warming and heat to 80°C) 配制储备液