


Product Data Sheet

Cas No.:	95-16-9	Cat. No:	PL06351
Product Name:	Benzothiazole		
Product synonym:	苯并噻唑;1,3-硫氮杂茛;间氮硫杂茛;苯硫氮醇;1,3-硫氮杂茛;间氮硫杂茛;Benzothiazole 苯并噻唑;苯并噻唑BT;中间体 苯并噻唑 BT		
Chemical name:	Benzothiazole		
MF:	C7H5NS	FW:	135.186300039291
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	N1C2C(=CC=CC=2)SC=1		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Benzothiazole 是一种天然的杂环核。Benzothiazole 杂环核具有抗癌、抗菌、抗糖尿病、抗炎、抗病毒、抗利什曼等多种生物活性。

生物活性	Benzothiazole is a natural occurring heterocyclic nuclei. Benzothiazole nucleus possesses a number of biological activities such as anticancer, antimicrobial, antidiabetic, anti-inflammatory, antileishmanial, and antiviral.
包装储存	Pure form -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 110 mg/mL (813.73 mM; Need ultrasonic)配制储备液