


Product Data Sheet

Cas No.:	80-65-9	Cat. No:	PL06357
Product Name:	3-Amino-2-oxazolidinone		
Product synonym:	3-氨基-2-恶唑烷酮;3-氨基-2-噁唑烷酮(AOZ);3-氨基-2-噁唑烷酮;3-氨基-2-恶唑烷酮(AOZ);3-氨基-2-恶唑烷酮标准品;3-氨基噁唑烷-2-酮;3-氨基甲基-2-恶唑烷酮;AOZ(硝基呋喃代谢产物) 标准品;WITEGANF005呋喃唑酮AOZ标准品;呋喃唑酮AOZ标准品;呋喃唑酮代谢物 标准品;硝基呋喃代谢产物-呋喃唑酮AOZ;3-氨-2-恶唑烷酮;3-氨基-2-恶唑啉酮;AOZ,呋喃唑酮代谢物 标准品;AOZ【呋喃唑酮代谢物,3-氨基-2-恶唑烷酮】		
Chemical name:	3-Amino-2-oxazolidinone		
MF:	C3H6N2O2	FW:	102.0919
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O1C(N(C([H]))([H])C1([H])([H])N([H])([H])=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

3-Amino-2-oxazolidinone (AOZ) 是呋喃唑酮 (HY-B1336) 的代谢物 (metabolite)。3-Amino-2-oxazolidinone 作为检测呋喃唑酮残留的一个指标。

生物活性	3-Amino-2-oxazolidinone (AOZ) is the metabolite of Furazolidone (HY-B1336). 3-Amino-2-oxazolidinone is always be detected as a indicator of furazolidone residues in vivo.
包装储存	4°C, protect from light In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)
溶解度数据	In Vitro: DMSO : 25 mg/mL (244.88 mM; ultrasonic and warming and heat to 60°C)配制储备液