


Product Data Sheet

Cas No.:	238742-88-6	Cat. No:	PL06360
Product Name:	Fmoc-D-Phe(4-CF3)-OH		
Product synonym:	Fmoc-D-4-三氟甲基苯丙氨酸;N-(9-芴甲氧羰基)-D-4-三氟甲基苯丙氨酸;5-氨基-1-叔丁基-1H-吡咯-3-甲腈;FMOC-4-(三氟甲基)-D-苯丙氨酸;Fmoc-4-(三氟甲基)-D-苯丙氨酸		
Chemical name:	Fmoc-D-Phe(4-CF3)-OH		
MF:	C25H20F3NO4	FW:	455.4258
Purity:	-	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	FC(C1C([H])=C([H])C=C([H])C=1[H])C([H])([H])[C@]([H])(C(=O)O[H])N([H])C(=O)OC([H])([H])C1([H])C2=C([H])C([H])=C([H])C([H])=C2C2=C([H])C([H])=C([H])C([H])=C12)(F)F		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Fmoc-D-Phe(4-CF₃)-OH 是苯丙氨酸衍生物。Fmoc-D-Phe(4-CF₃)-OH 可用于蛋白-蛋白相互作用的肽抑制剂的研宄。

生物活性	Fmoc-D-Phe(4-CF ₃)-OH is Phenylalanine derivative. Fmoc-D-Phe(4-CF ₃)-OH can be used for the research of peptide inhibitors of protein-protein interactions.
包装储存	Powder -20°C 3 years; 4°C 2 years