


Product Data Sheet

Cas No.:	166247-63-8	Cat. No:	PL06509
Product Name:	7-Deaza-2'-deoxy-7-iodoadenosine		
Product synonym:	7-脱氮-2'-脱氧-7-碘腺苷;7-脱氮-2'-脱氧-7- IODOADENOSINE;7-Deaza-7-I-2'-dA		
Chemical name:	7-Deaza-2'-deoxy-7-iodoadenosine		
MF:	C11H13IN4O3	FW:	376.1504
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	IC1C2=C(N([H])[H])N=C([H])N=C2N(C=1[H])[C@@]1([H])C([H])([H])[C@@]([H])([C@@]([H])(C([H])([H])O[H])O1)O[H]		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Deaza-2'-deoxy-7-iodoadenosine 是包含 7-Deazaadenine 的修饰的寡核苷酸。

生物活性	7-Deaza-2-deoxy-7-iodoadenosine is a modified oligonucleotide containing 7-Deazaadenine.
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 50 mg/mL (132.93 mM; Need ultrasonic)配制储备液