


Product Data Sheet

Cas No.:	1869126-64-6	Cat. No:	PL06562
Product Name:	Azido-PEG4-Val-Cit-PAB-MMAE		
Product synonym:	-		
Chemical name:	Azido-PEG4-Val-Cit-PAB-MMAE		
MF:	C69H113N13O17	FW:	1396.71
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	[C@@H]([C@]1([H])CCCN1C(=O)C[C@@H](OC)[C@H]([C@@H](C)CC)N(C)C(=O)[C@H](C(C)C)NC(=O)[C@H](C(C)C)N(C)C(=O)OCC1C=CC(=CC=1)NC(=O)[C@H](CCCN(C(=O)N)NC(=O)[C@H](C(C)C)NC(=O)CCOCCOCCOCCOCCN=[N+]=[N-])(OC)[C@@H](C)C(=O)N[C@H](C)[C@H](C1C=CC=CC=1)O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Azido-PEG4-Val-Cit-PAB-MMAE 是抗体偶联药物的一部分 (drug-linker conjugate for ADC)。由多肽 linker 和抗有丝分裂剂 MMAE (一种 tubulin 抑制剂) 连接而成。

生物活性	Azido-PEG4-Val-Cit-PAB-MMAE is a drug-linker conjugate for ADC by using the anti-mitotic agent, monomethyl auristatin E (MMAE, a tubulin inhibitor), linked via the cleavable linker Azido-PEG4-Val-Cit-PAB-OH.
IC50 & Target[1][2]	Drug-Linker Conjugates for ADC
包装储存	-20°C, protect from light, stored under nitrogen 该产品在溶液状态不稳定, 建议您现用现配, 即刻使用。