


Product Data Sheet

Cas No.:	1547162-46-8	Cat. No:	PL06610
Product Name:	CC-17369		
Product synonym:	-		
Chemical name:	CC-17369		
MF:	C13H11N3O5	FW:	289.243542909622
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C1C(CCC(N1)=O)N1C(C2C(=CC=C(C=2C1=O)N)O)=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

CC-17369 (7-Hydroxy pomalidomide) 是 Pomalidomide 的代谢物。CC-17369 基于泊马利度胺的cereblon (CRBN) 配体。CC-17369 可通过连接物连接到蛋白质配体合成 PROTAC。

生物活性	CC-17369 (7-Hydroxy pomalidomide) is a metabolite of Pomalidomide. CC-17369 is the Pomalidomide -based cereblon (CRBN) ligand used in the recruitment of CRBN protein. CC-17369 can be connected to the ligand for protein by a linker to form PROTAC.
IC50 & Target[1][2]	Cereblon
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 50 mg/mL (172.87 mM; Need ultrasonic)配制储备液