


Product Data Sheet

Cas No.:	1570374-32-1	Cat. No:	PL06937
Product Name:	TAK1-IN-4		
Product synonym:	-		
Chemical name:	TAK1-IN-4		
MF:	C18H17N3O3	FW:	323.345884084702
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C(C1C=CC=C(C(=O)O)C=1)NC1=NC2C=CC=CC=2N1CCC		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

TAK1-IN-4 (Compound 14) 是一种 TAK1 抑制剂。

生物活性	TAK1-IN-4 (Compound 14) is a TAK1 inhibitor.
IC50 & Target[1][2]	TAK1
包装储存	4°C, protect from light In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)
溶解度数据	In Vitro: DMSO : 83.33 mg/mL (257.71 mM; Need ultrasonic)配制储备液