

体外研究(In Vitro)	o-Phenanthroline (1,10-Phenanthroline, 1mM) inhibits autolysis of sea cucumber body wall. has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	4°C, sealed storage, away from moisture In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)
溶解度数据	In Vitro: DMSO : 100 mg/mL (504.49 mM; Need ultrasonic)H ₂ O : 2 mg/mL (10.09 mM; ultrasonic and warming and heat to 60°C)配制储备液