


Product Data Sheet

Cas No.:	92434-13-4	Cat. No:	PL07611
Product Name:	NHE3-IN-2		
Product synonym:	-		
Chemical name:	NHE3-IN-2		
MF:	C15H12N5CL	FW:	297.742
Purity:	≥95%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	NC(NC1=NC(C2=CC=CC=C2)=C3C=C(Cl)C=CC3=N1)=N		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

NHE3-IN-2 是一种 Na⁺/H⁺ 交换体 3 (NHE3) 抑制剂 (详细信息请参考专利 WO2001079186A1 6-Chlor-4-phenyl-2-chinazolinyl-guanidin)。

生物活性	NHE3-IN-2 is a Na/H exchanger-3 (NHE3) inhibitor (patent WO2001079186A1, example 6-Chlor-4-phenyl-2-chinazolinyl-guanidin).
包装储存	4°C, protect from light In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)
溶解度数据	In Vitro: DMSO : 3.33 mg/mL (11.18 mM; ultrasonic and warming and heat to 60°C)配制储备液