


Product Data Sheet

Cas No.:	876707-99-2	Cat. No:	PL07973
Product Name:	Nucleoside-Analog-1		
Product synonym:	-		
Chemical name:	Nucleoside-Analog-1		
MF:	C9H9N5O5	FW:	267.198
Purity:	≥95%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O[C@@H]1[C@](CO)(N=[N+]=[N-])O[C@@H](N2C=C3)[C@H]1OC2=NC3=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Nucleoside-Analog-1是抗病毒核苷R1479的类似物，用于HCV 抗病毒核苷类似物的合成。

生物活性	Nucleoside-Analog-1 is a 4'-Azidocytidine analogue against Hepatitis C virus replication.
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: H ₂ O : 1.89 mg/mL (7.07 mM; Need ultrasonic)配制储备液