


## Product Data Sheet

Cas No.:	229002-10-2	Cat. No:	PL08741
Product Name:	MAPK13-IN-1		
Product synonym:	-		
Chemical name:	MAPK13-IN-1		
MF:	C20H23N5O2	FW:	365.42892
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C(NC1=CC=C(OC2=CC=NC=C2)C=C1)NC3=CC(C(C)(C)C)=NN3C		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

## Product Description

MPAK13-IN-1 是 MAPK13 (p38δ) 的一个抑制剂, 其 IC50 值为 620 nM。

生物活性	MPAK13-IN-1 is a MAPK13 (p38δ) inhibitor, with an IC 50 of 620 nM.
IC50 & Target[1][2]	MAPK13 (p38δ) 620 nM (IC50)
体外研究(In Vitro)	MPAK13-IN-1 exhibits an IC50 of 4.63 μM in vero E6 cells. has not independently confirmed the accuracy of these methods. They are for reference only.Cell Viability Assay
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : ≥ 250 mg/mL (684.13 mM)配制储备液