

## Product Data Sheet

Cas No.:		Cat. No:	PL08978
Product Name:	Human PD-L1 inhibitor V TFA		
Product synonym:	-		
Chemical name:	Human PD-L1 inhibitor V TFA		
MF:	-	FW:	-
Purity:	≥96%	Batch No.:	-
Storage:			
λmax:	-	Formulation:	-
Solubility :	-		
SMILES :	-		
InChI Code:	-		
InChI Key:	-		
WARNING This product is not for human or veterinary use.			

## Product Description

Human PD-L1 inhibitor V TFA 是一种 human PD-1 protein 结合肽，K<sub>d</sub> 值为 3.32 μM，可抑制人 PD-1/PD-L1 的结合。

生物活性	Human PD-L1 inhibitor V TFA, a human PD-1 protein binding peptide with a K <sub>d</sub> value of 3.32 μM. Human PD-L1 inhibitor V TFA inhibit the interaction of hPD-1/hPD-L1.
体外研究(In Vitro)	Human PD-L1 inhibitor V (Ar5M_2) is a peptide designed with anchor residues Y56, R113, A121, D122 and Y123. has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	Sealed storage, away from moisturePowder -80°C 2 years; -20°C 1 year
溶解度数据	In Vitro: H <sub>2</sub> O : 100 mg/mL (62.51 mM; Need ultrasonic)配制储备液