

Product Data Sheet

Cas No.:		Cat. No:	PL09866
Product Name:	PKI 14-22 amide,myristoylated TFA		
Product synonym:	-		
Chemical name:	PKI 14-22 amide,myristoylated TFA		
MF:	-	FW:	-
Purity:	≥99%	Batch No.:	-
Storage:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	-		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

PKI 14-22 amide,myristoylated TFA 是一种有效的 cAMP 依赖性 PKA 抑制剂。PKI 14-22 amide,myristoylated TFA 可降低 IgG 介导的吞噬反应，也可抑制中性粒细胞粘附。

生物活性	PKI 14-22 amide,myristoylated TFA is a potent cAMP-dependent PKA inhibitor. PKI 14-22 amide,myristoylated TFA reduces the IgG-mediated phagocytic response and also inhibits neutrophil adhesion.
包装储存	Sealed storage, away from moisture Powder -80°C 2 years; -20°C 1 year
溶解度数据	In Vitro: H ₂ O : 1.85 mg/mL (1.40 mM; Need ultrasonic)配制储备液