

Product Data Sheet

Cas No.:		Cat. No:	PL11246
Product Name:	Dapoxetine-d7 (hydrochloride)		
Product synonym:	-		
Chemical name:	Dapoxetine-d7 (hydrochloride)		
MF:	-	FW:	-
Purity:	≥99%	Batch No.:	-
Storage:			
λmax:	-	Formulation:	-
Solubility :	-		
SMILES :	-		
InChI Code:	-		
InChI Key:	-		
WARNING This product is not for human or veterinary use.			

Product Description

Dapoxetine-d₇ (hydrochloride) 是 Dapoxetine hydrochloride 的氘代物。Dapoxetine hydrochloride 是一种短效的选择性 5-羟色胺再摄取抑制剂 (SSRI)。

生物活性	Dapoxetine-d 7 (hydrochloride) is the deuterium labeled Dapoxetine hydrochloride. Dapoxetine hydrochloride is a short-acting selective serotonin reuptake inhibitor (SSRI)[1].
包装储存	4°C, sealed storage, away from moistur In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)
溶解度数据	In Vitro: DMSO : 250 mg/mL (716.50 mM; Need ultrasonic)H ₂ O : 50 mg/mL (143.30 mM; Need ultrasonic)配制储备液