

Product Data Sheet

Cas No.:		Cat. No:	PL11573
Product Name:	Acetyl-PHF6 amide TFA		
Product synonym:	-		
Chemical name:	Acetyl-PHF6 amide TFA		
MF:	-	FW:	-
Purity:	≥95%	Batch No.:	-
Storage:			
λmax:	-	Formulation:	-
Solubility :	-		
SMILES :	-		
InChI Code:	-		
InChI Key:	-		
WARNING This product is not for human or veterinary use.			

Product Description

Acetyl-PHF6 amide TFA (AcPHF6 TFA) 是 tau 衍生的六肽。

生物活性	Acetyl-PHF6 amide TFA (AcPHF6 TFA) is a tau derived hexapeptide.
包装储存	Sealed storage, away from moisture and light Powder -80°C 2 years; -20°C 1 year
溶解度数据	In Vitro: DMSO : 67.5 mg/mL (74.67 mM; Need ultrasonic) 配制储备液