


Product Data Sheet

Cas No.:	163018-26-6	Cat. No:	PL12545
Product Name:	3-O-(2-Aminoethyl)-25-hydroxyvitamin D3		
Product synonym:	25 - 羟基维生素D33,3' - 氨基丙基醚;25 - 羟基维生素D33,3'' - 氨基丙基醚;25-羟基维生素 d3 3,3;5-羟基-1H-吡啶-2-酮;聚(苯乙烯-马来酸酐)部分末端丙酯,异丙苯;3-O-(2-氨基丙基)-25-羟基维他命 D3;25 - 羟基维生素D33		
Chemical name:	3-O-(2-Aminoethyl)-25-hydroxyvitamin D3		
MF:	C30H51NO2	FW:	457.7314
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	<chem>O([H])C(C([H])([H]))(C([H])([H])[H])C([H])([H])C([H])([H])C([H])([H])C([H])([H])C@ @]([H])(C([H])([H])[H])C@ @]1([H])C([H])([H])C([H])([H])C@ @]2([H])/C=C\([H])/C/[H])=C3C(=C([H])([H])C([H])([H])C([H])([H])C@ @]([H])(C3([H])[H])OC([H])([H])C([H])([H])C([H])([H])N([H])([H])/C([H])([H])C([H])([H])C([H])([H])C@ @]12C([H])([H])[H]</chem>		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

3-O-(2-Aminoethyl)-25-hydroxyvitamin D3 是一种Vitamin D3 衍生物。

生物活性	3-O-(2-Aminoethyl)-25-hydroxyvitamin D3 is a Vitamin D3 derivative.
IC50 & Target[1][2]	Vitamin D/Vitamin D receptor
体外研究(In Vitro)	3-O-(2-Aminoethyl)-25-hydroxyvitamin D3 is a 25-hydroxyvitamin D3 derivative. has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	4°C, protect from light, stored under nitrogen 该产品在溶液状态不稳定，建议您现用现配，即刻使用。
溶解度数据	In Vitro: DMSO : ≥ 100 mg/mL (218.47 mM)配制储备液