


Product Data Sheet

Cas No.:	1641940-94-4	Cat. No:	PL12540
Product Name:	Alfacalcidol-d6		
Product synonym:	阿尔法骨化醇-D6		
Chemical name:	Alfacalcidol-d6		
MF:	C27H44O2	FW:	406.6740
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O([H])[C@@]1([H])C([H])([H])[C@@]([H])(C=C([H])[H])/C/C1([H])[H]=C/[H]\C([H])=C1/C([H])([H])C([H])([H])C([H])([H])[C@@]2(C([H])([H])[H])[C@@]([H])C([H])([H])C([H])([H])[C@]2([H])[C@]([H])(C([H])([H])[H]C([H])([H])C([H])([H])C([H])([H])C([H])([2H])([2H])C([2H])([2H])[2H])O[H]		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Alfacalcidol-d₆ 是 Alfacalcidol 的氘代物。Alfacalcidol 是非选择性的 VDR 活化剂。

生物活性	Alfacalcidol-d 6 is the deuterium labeled Alfacalcidol. Alfacalcidol is a non-selective VDR activator.
包装储存	Powder -20°C 3 years; 4°C 2 years