


Product Data Sheet

Cas No.:	118584-54-6	Cat. No:	PL12556
Product Name:	VD3-d6		
Product synonym:	维生素 D3-D6;VD3-D6		
Chemical name:	VD3-d6		
MF:	C27H44O	FW:	384.63766
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O[C@@H](CCC/C=C)CC1=C\C=C2[C@@H](CC[C@@H]3[C@@H](CCCC(C(2H))([2H])[2H])C([2H])([2H])[2H])C([H])([C@]3(C)CCC\2		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

VD3-d₆ 是维生素 D3 的氘代物，可作为人血清中维生素 D3 代谢物测定工具。

生物活性	VD3-d 6 is the deuterium labeled Vitamin D3; tools for determination of Vitamin D3 metabolites in human serum.
包装储存	-20°C, protect from light, stored under nitrogen 该产品在溶液状态不稳定，建议您现用现配，即刻使用。
溶解度数据	In Vitro: DMSO : 100 mg/mL (255.97 mM; Need ultrasonic)配制储备液