

Product Data Sheet

[illegible]

Product Description

Tween 80 (Polysorbate 80) 是一种表面活性剂，已被广泛用作药理学实验的溶剂。Tween 80 可以减少细菌附着并抑制生物膜的形成。

生物活性	Tween 80 (Polysorbate 80), a surfactant, has been widely used as a solvent for pharmacological experiments. Tween 80 can also reduce bacterial attachment and inhibit biofilm formation.
------	--

体外研究(In Vitro)	Tween 80 (Polysorbate 80) is a surfactant which can also reduce bacterial attachment and inhibit biofilm formation. Solubility enhancement as high as 100-fold is achieved using 1% Tween 80 in vitro . has not independently confirmed the accuracy of these methods. They are for reference only.
体内研究(In Vivo)	Tween 80 (Polysorbate 80) (IP; 10 ml/kg) significantly decreases locomotor activity at a concentration of 32% in saline (v/v) in male outbred (BALB/c x DBA/Z)F1 mice weighing 24-30 g. has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	Pure form -20°C 3 years; 4°C 2 years