


Product Data Sheet

Cas No.:	2380228-45-3	Cat. No:	PL13220
Product Name:	DC-TEADin02		
Product synonym:	-		
Chemical name:	DC-TEADin02		
MF:	C19H17NO2S	FW:	323.41
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	C(S(NC1=CC=C(CC2=CC=C3C(=C2)C=CC=C3)C=C1)(=O)=O)=C		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

DC-TEADin02 是一种有效的 TEAD 自棕榈酰化抑制剂。DC-TEADin02 具有 TEAD 自棕榈酰化抑制作用，IC₅₀ 值为 197 nM。DC-TEADin02 可用于发育、再生和组织稳态研究。

生物活性	DC-TEADin02 is a potent TEAD autopalmitylation inhibitor. DC-TEADin02 has TEAD autopalmitylation inhibitory with the IC ₅₀ value of 197 nM. DC-TEADin02 can be used for the research of development, regeneration and tissue homeostasis.
体外研究(In Vitro)	DC-TEADin02 具有 TEAD 自棕榈酰化抑制作用，IC ₅₀ 值为 197 nM。 DC-TEADin02 (0-25 μM) 抑制 TEAD 转录活性，导致 YAP 相关下游基因表达下调。 has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	4°C, protect from light In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)
溶解度数据	In Vitro: DMSO : 100 mg/mL (309.21 mM; Need ultrasonic)配制储备液