


Product Data Sheet

Cas No.:	2748749-29-1	Cat. No:	PL13357
Product Name:	IL-17A modulator-1		
Product synonym:	-		
Chemical name:	IL-17A modulator-1		
MF:	C33H31N5O4	FW:	561.63
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	C(C1C=CC=CC=1)(C1C=CC=CC=1)[C@@H](C(=O)NC1C=CC(C2=C(C(=NC=C2OC)OC)=CC=1)NC(C1N(N=CC=1)C)=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

IL-17A modulator-1 是一种 IL-17A 调节剂，来自专利 WO2021239743+A1，example 9。IL-17A modulator-1 抑制 IL-17A 的生物学效应，pIC₅₀ 为 8.2。IL-17A modulator-1 可用于研究与调节 IL-17A 活性相关的疾病或病症，包括具有免疫成分或自身免疫病理的疾病、癌症和神经退行性疾病。

生物活性	IL-17A modulator-1 is a IL-17A modulator, extracted from patent WO2021239743+A1, example 9. IL-17A modulator-1 inhibits the biological action of IL-17A with a pIC ₅₀ of 8.2. IL-17A modulator-1 can be used for the research of diseases or disorders associated with modulation of IL-17A activity including diseases with an immune component or autoimmune pathol, cancer and neurodegenerative disorders.
IC50 & Target[1][2]	IL-17A 8.2 (pIC ₅₀)
包装储存	4°C, sealed storage, away from moisture and light In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)
溶解度数据	In Vitro: DMSO : 100 mg/mL (178.05 mM; Need ultrasonic)配制储备液