


Product Data Sheet

Cas No.:	2468048-24-8	Cat. No:	PL13441
Product Name:	FHT-1204		
Product synonym:	-		
Chemical name:	FHT-1204		
MF:	C24H23N5O5S2	FW:	525.60
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	N1(S(C)(=O)=O)C=CC(N[C@@H](COC)C(=O)NC2=NC(C3=CC=CC(C4C=CN=CC=4)=C3)=CS2)=O)=C1		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

FHT-1204 是一种有效的 SMARCA4/SMARCA2 ATPase (BRG1 和 BRM) 抑制剂, IC₅₀ ≤10 nM (WO2020160180A1; 化合物 70)。

生物活性	FHT-1204 is a potent SMARCA4/SMARCA2 ATPase (BRG1 and BRM) inhibitor with IC 50 s of ≤10 nM (WO2020160180A1; compound 70).
包装储存	Powder -20°C 3 years; 4°C 2 years
溶解度数据	In Vitro: DMSO : 25 mg/mL (47.56 mM; ultrasonic and warming and heat to 60°C)配制储备液