


Product Data Sheet

Cas No.:	96801-55-7	Cat. No:	PL14005
Product Name:	2-Thio-PAF		
Product synonym:	1-O-十六烷基-2-硫代乙酰基-2-脱氧-SN-甘油基-3-磷酸胆碱		
Chemical name:	2-Thio-PAF		
MF:	C26H55NO6PS+	FW:	540.7568
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	CCCCCCCCCCCCCCCCOC[C@H](COP(=O)([O-])OCC[N+](C)(C)C)SC(=O)C		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

2-Thio-PAF 是 PAF 的合成类似物，用于 PAF 乙酰水解酶的比色测定。

生物活性	2-Thio-PAF, a synthetic analog of PAF, is used in a colorimetric assay for PAF acetylhydrolases.
体外研究(In Vitro)	2-Thio-PAF induced platelet-aggregation that was inhibited by previous incubation of the substrate with salivary gland homogenates or saliva. has not independently confirmed the accuracy of these methods. They are for reference only.
包装储存	Solution, -20°C, 2 years