


Product Data Sheet

Cas No.:	349554-02-5	Cat. No:	PL14775
Product Name:	Enalapril-d5 (maleate)		
Product synonym:	马来酸依那普利-D5;马来酸依那普利-D5氘代内标;依那普利D5;依那普利-d5马来酸盐(苯基-d5);依那普利,马来酸D5;马来酸依那普利D5		
Chemical name:	Enalapril-d5 (maleate)		
MF:	C20H23D5N2O5.C4H4O4	FW:	487.47918
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C(O)[C@H]1N(C([C@H](C)N[C@H](C(OCC)=O)CCC2=C([2H])C([2H])=C([2H])C([2H])=C2[2H])=O)CCC1.O=C(O)/C=C/C(O)=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Enalapril-d₅ (maleate) 是 Enalapril 的氘代物。

生物活性	Enalapril-d 5 (maleate) is deuterium labeled Enalapril, which is an angiotensin converting enzyme (ACE) inhibitor.
包装储存	-20°C, sealed storage, away from moistur In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)