


Product Data Sheet

Cas No.:	349554-00-3	Cat. No:	PL14791
Product Name:	Enalaprilat-d5		
Product synonym:	依那普利-D5;依那普利拉-D5;依那普利拉-d5(苯基-d5);依那普利拉-D5氘代内标;依那普利拉 D5;依那普利拉特D5;[2H5]-ENALAPRILAT;ENALAPRILAT: H2O		
Chemical name:	Enalaprilat-d5		
MF:	C18H24N2O5	FW:	353.424373626709
Purity:	≥99%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	O=C([C@]([2H])(C([2H])([2H])N([2H])(C@H)(C(=O)O)CCC1C=CC=CC=1)N1CCC[C@@]1([2H])C(=O)O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

Enalaprilat-d₅ 是 Enalaprilat 的氘代物。

生物活性	Enalaprilat-d 5 is the deuterium labeled Enalaprilat(MK-422), which is an angiotensin-converting enzyme (ACE) inhibitor.
包装储存	Powder -20°C 3 years; 4°C 2 years