


Product Data Sheet

Cas No.:	62288-83-9	Cat. No:	PC10603
Product Name:	Desmopressin Acetate		
Product synonym:	1-去氨基-8-右旋精氨酸加压素;醋酸去氨加压素;去氨加压素		
Chemical name:	Desmopressin Acetate		
MF:	C48H68N14O14S2	FW:	1129.2689
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	S1C([H])([H])[C@@]([H])(C(N2C([H])([H])C([H])([H])C([H])([H])[C@]2([H])C(N([H])[C@]([H])(C(N([H])C([H])([H])C(N([H])([H])=O)=O)C([H])([H])C([H])([H])C([H])([H])/N=C\N([H])([H])/N([H])([H])=O)=O)N([H])C([C@]([H])([H])(C([H])([H])C(N([H])([H])=O)N([H])C([C@]([H])([H])(C([H])([H])C(N([H])([H])=O)N([H])C([C@]([H])([H])(C2C([H])=C([H])C([H])=C([H])C=2[H])N([H])C([C@]([H])([H])([H])C2C([H])=C([H])C=C([H])C=2[H])O[H])N([H])C(C([H])([H])C([H])([H])S1)=O)=O)=O)=O.O([H])C(C([H])([H])([H])=O		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

精氨酸加压素的合成类似物,Desmopressin(DDAVP)醋酸盐是抗利尿激素精氨酸加压素类似物。