


Product Data Sheet

Cas No.:	2480-28-6	Cat. No:	PC12000
Product Name:	N2-Methyl-L-arginine		
Product synonym:	-		
Chemical name:	N2-Methyl-L-arginine		
MF:	C7H16N4O2	FW:	188.22754
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:			
λmax:	-	Formulation:	-
Solubility :			
SMILES :	CNC(C(O)=O)CCCNC(=N)N		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

L-精氨酸摄取抑制剂