


Product Data Sheet

Cas No.:	10347-81-6	Cat. No:	PC10703
Product Name:	Maprotiline HCl		
Product synonym:	盐酸马普替林;马普替林;吗丙啉;路滴美;马替普林盐酸盐;Maprotiline Hydrochloride 马替普林盐酸盐;黑染料;马普替林盐酸;麦普替林-D3盐酸;盐酸马普替林 EP标准品;盐酸马普替林 USP标准品;盐酸马普替林 标准品;1-(3-甲基氨基丙基)二苯并[b,e]二环[2.2.2]辛二烯盐酸盐;9-(γ-氨基丙基)-9,10-二氢-9,10-乙醇蒽盐酸盐;盐酸马替普林		
Chemical name:	Maprotiline HCl		
MF:	C20H24CLN	FW:	313.8643
Purity:	≥98%	Batch No.:	-
Storage:			
Structural formula:	 HCl		
λmax:	-	Formulation:	-
Solubility :			
SMILES :	Cl[H].N([H])(C([H])([H])([H])C([H])([H])C([H])([H])C([H])([H])C12C3=C([H])C([H])=C([H])C([H])=C3C([H])(C3=C([H])C([H])=C([H])C([H])=C13)C([H])([H])C2([H])[H])		
InChI Code:	-		
InChI Key:			
WARNING This product is not for human or veterinary use.			

Product Description

选择性去甲肾上腺素再摄取抑制剂,Maprotiline盐酸盐是去甲肾上腺素重吸收抑制剂, 为四环抗抑郁化合物 (TeCA) 。